

## LABNOTE

### ELECTRON-DENSITY DETERMINATION OF A BOROXINE COMPOUND FROM HIGH-RESOLUTION X-RAY DATA MEASURED WITH A STOE IPDS 2T


1 | 2

Simon Grabowsky<sup>1</sup>, Stephan D. Hoffmann<sup>2</sup>, Usman Ahmad<sup>1</sup>, Jens Beckmann<sup>1</sup>, Jens Richter<sup>2</sup> and Peter Luger<sup>1</sup>.  
<sup>1</sup>Freie Universität Berlin, Institut für Chemie und Biochemie / Anorganische Chemie, Fabeckstr. 36a, 14195 Berlin.  
<sup>2</sup>Stoe & Cie GmbH, Hilpertstr. 10, 64295 Darmstadt.


The boroxine compound shown in Figure 1 originally served as a model compound to exemplify concepts of ringstrain and electron deficiency.<sup>[1]</sup>

Its electron-density (ED) distribution was derived from a high-resolution X-ray data set of 286975 reflections collected at 90 K on an IPDS 2T that allowed  $2\theta$ -movement to reach high regions in reciprocal space. A typical frame of the measurement is presented in Figure 2.

3082 frames with an exposure time of one minute were collected, using  $\omega$ -scans with an  $\omega$  increment of 0.3 degrees. To gain a better signal to noise ratio for the weak reflections at high  $2\theta$  angles, a second data set of 1932 frames with an exposure time of five minutes was collected. Both data sets were scaled together and used for the ED determination.



**Fig. 1:** Molecular structure and electrostatic potential (esp) of the examined boroxine compound.



**Fig. 2:** Typical frame of the measurement with an IPDS 2T


Hence, a maximum resolution of  $\sin\theta_{\max}/\lambda = 1.22 \text{ \AA}^{-1}$  ( $2\theta_{\max} = 120^\circ$ ), an overall completeness of 97.5% and a redundancy of 9.1 could be reached. The data set, reduced to 30876 unique reflections ( $R_{\text{int}} = 4.1\%$ ), was properly suited to model successfully the ED distribution using the multipole formalism as implemented in XD2006<sup>[2]</sup>: 778 parameters,  $wR_{\text{mult}}(F^2) = 1.5\%$ ,  $\text{GoF} = 2.8$ .

## LABNOTE


### ELECTRON-DENSITY DETERMINATION OF A BOROXINE COMPOUND FROM HIGH-RESOLUTION X-RAY DATA MEASURED WITH A STOE IPDS 2T

2 | 2

Figures 1, 3 and 4 show the esp and deformation densities as first ED results. An intra-molecular hydrogen bond C11-H11...O3 could be found that exhibits a strong electrostatic nature, see Fig. 1. The deformation densities show the shape of bonding and non-bonding effects. For example, the oxygen lone pairs are clearly developed in the boroxine plane, see Fig. 4.



**Fig. 3:** Static deformation density map in the plane of the pyridine ring.



**Fig. 4:** Static deformation density map in the plane of the boroxine ring.

As conclusion, the STOE IPDS 2T has proven as a powerful instrument in the collection of high resolution and high quality data for electron density determination.

[1] J. Beckmann, D. Dakternieks, A. Duthie, A. E. K. Lim, E. R. T. Tiekink, *J. Organomet. Chem.* 2001, 633, 149.

[2] A. Volkov, P. Macchi, L. J. Farrugia, C. Gatti, P. Mallinson, T. Richter, T. Koritsanszky, XD2006