

LABNOTE IPDS 2T CHEMICAL CRYSTALLOGRAPHY APPLICATIONS WITH THE STOE IPDS 2T DIFFRACTOMETER EQUIPPED WITH A XENOCs GENIX CU K_{α} MICROBEAM DELIVERY SYSTEM

*K. Schenk, C. Besnard, M. Schiltz. Laboratoire de Cristallographie.
École Polytechnique Fédérale de Lausanne (EPFL). Switzerland.*

CHARACTERISTICS




Fig.1: IPDS2T-GeniX setup at the EPFL

Although Mo K_{α} radiation is widely used in routine small molecule structure determination, a number of applications in chemical crystallography are best treated with Cu K_{α} radiation. These include studies on weakly diffracting and/or very small crystals, measurements on partially disordered or fibrous samples, absolute structure determination and others.

To evaluate the suitability of the STOE image-plate system IPDS2T equipped with a GeniX Cu K_{α} microbeam system manufactured by XENOCs for applications in chemical crystallography, we have carried out trial measurements on two samples that are commonly used for testing purposes.

1. Data collection on YLID (C₁₁H₁₀O₂S)

YLID is a small-molecule standard widely used in the single-crystal diffraction. It belongs to the orthorhombic system so that the refinement of the cell angles provides a good check of the alignment, geometric accuracy and general performances of a laboratory diffractometer setup.

The GeniX microbeam system was operated at 50 kV, 1 mA. The circular image plate (scanned surface: 170 mm radius) was positioned at a distance of 40 mm from the sample with a 60° 2 θ offset, thus allowing data to be collected to minimum Bragg spacings of 0.83 Å. A highly redundant data set, consisting of three 0-180° omega-scans recorded at three different phi angles (0°, 45° and 90°) was collected. The data comprise a total of 540 frames, each corresponding to a 1° omega rotation and an exposure of 60 seconds. The frames were indexed and integrated with the X-Area software from STOE. Data reduction (scaling and merging) was carried out with the program XRed from STOE. For comparison, reference data on the same sample were also collected on a sealed tube Mo K_{α} source in our laboratory, using an IPDS 2 diffractometer from STOE. Structure refinement was carried out with SHELXL.


Refined cell parameters are reported in Table 1 for both measurements, along with available reference data (G.L.Bryant¹). Data reduction and structure refinement statistics are reported in Table 2. The data collected on the STOE IPDS 2T / Xenocs-GeniX system are slightly better than the literature values that were obtained from one of the best reputed Cu K_{α} instruments on the market (Bryant¹). Most certainly, they would be even better had we been given a better mounted and shaped test-crystal. The comparison of the two measurements ascertains the superior intensity of the Xenocs/GeniX source. For roughly the same R σ value the exposure time with Mo radiation (53kV, 43mA) had to be five times bigger than that with the microbeam system.

¹ <http://www.avtechlabs.com/ylid2/ylid.htm>

LABNOTE IPDS 2T CHEMICAL CRYSTALLOGRAPHY APPLICATIONS WITH THE STOE IPDS 2T DIFFRACTOMETER EQUIPPED WITH A XENOCs GENIX CU K_α MICROBEAM DELIVERY SYSTEM

2 | 3

2. Data Collection on Zeolite LTA $\{Na_{12}[Al_{12}Si_{12}O_{48}] \cdot 27 H_2O\}_8$



Zeolites are aluminosilicate minerals. Their porous molecular structures have made them very attractive materials for a wide range of industrial applications including molecular sieves, catalysts and detergents.

Zeolite A exhibits the LTA (Linde Type A) structure. The *topological* unit cell of symmetry is doubled by the chemical partition between Al and Si, giving rise to a *supercell* with unit-cell parameters of 24.572 Å. The doubling of the unit-cell generates very weak superlattice reflections.

The STOE image plate system *IPDS 2T*, equipped with a *GeniX* microbeam system from *XENOCs* was used to collect diffraction data on a small crystal of zeolite A.

Fig.2: The framework structure of Zeolite A with Si-O and Al-O bonds colored respectively in blue and red

The GeniX microbeam system was operated at 50kV, 1mA and an exposure time of 90 seconds per 1° omega-rotation allowed to record data with good statistics (*Fig. 3*). The image plate was positioned at a distance of 70 mm from the sample and was offset in 2θ by 45°. Three omega scans were collected at different phi angles producing highly redundant data. Data collection parameters and data reduction statistics are reported in Tables 3 and 4. In a reciprocal lattice reconstruction from the experimental data (*Fig. 4*), the superlattice reflections are clearly visible, thus substantiating the quality of the data. The framework structure (*Fig. 2*) could be solved using SHELXS with default parameters.




Fig.3: Diffraction frame recorded on a crystal of Zeolite A. Scan width 1° (omega), exposure time 90 sec. Crystal-detector distance 70 mm, 2θ offset 45°




Fig.4: Reciprocal lattice reconstruction of the (hk0) layer from experimental data. The weak superlattice reflections are clearly visible.

CONCLUSION:

Accurate data were collected on samples from organic molecules and complicated inorganic structures with the combined GeniX - IPDS 2T setup installed at EPFL-Lausanne.



LABNOTE IPDS 2T CHEMICAL CRYSTALLOGRAPHY APPLICATIONS WITH THE STOE IPDS 2T DIFFRACTOMETER EQUIPPED WITH A XENOC S GENIX CU K α MICROBEAM DELIVERY SYSTEM

3 | 3

Data		a [Å]	b [Å]	c [Å]	α [Å]	β [Å]	γ [Å]
IPDS 2T + GeniX Cu K α	nc	5.9603(5)	9.0354(5)	18.3859(12)	90.036(5)	90.033(6)	90.036(6)
	c	5.9677(4)	9.0336(5)	18.3883(11)			
IPDS 2 + sealed tube Mo K α	nc	5.9606(5)	9.0337(7)	18.3818(15)	90.014(7)	90.010(6)	90.027(6)
	c	5.9606(3)	9.0338(5)	18.3796(15)			
Reference (Bryant) Cu K α	c	5.96260(10)	9.03940(10)	18.3890(2)			

Table 1: Cell parameter refinement of YLID. nc: non constrained refinement, c: constrained refinement (i.e. cell angles are fixed at 90°)

Data	d_{\min} / Å	R_{int}	R_{σ}	SAV	IEq	R_1	Flack par.
IPDS2T + GeniX Cu K α	0.8354	0.0629	0.0362	0	291	0.0292	-0.019(21)
	0.9020	0.0627	0.0353	0	285	0.0240	-0.011(18)
IPDS2 + sealed tube Mo K α	0.5717	0.0255	0.0242	0	3	0.0363	-0.005(61)
Reference (Bryant et al.) Cu K α	0.9020	0.0644	0.0343	13	218	0.0277	-0.022(19)

Table 2: Data reduction and structure refinement statistics for Ylid. SAV: Systematic absence violations, IEq: Inconsistent Equivalents

Refined unit cell parameter: a / Å	24.5719 (0.0009)
Space group	$F m - 3 c$
Number of measured reflections	16249
Number of unique reflections	881
R_{int}	0.0575
Absorption correction transmission factors	Face-indexed 0.41 0.51
R_{int} after absorption correction	0.0545

Table 3: Data collection and reduction statistics for Zeolite A

d_{\min} / Å	14.19	2.60	1.90	1.62	1.42	1.31	1.21	1.13	1.06	0.99	0.93
Mean $I/\sigma(I)$	28.12	12.48	16.43	18.12	15.40	15.01	7.11	5.34	5.76	3.25	
R_{int}	28.12	0.076	0.068	0.054	0.043	0.040	0.050	0.054	0.072	0.152	
Completeness / %	98	100	100	100	100	100	100	100	100	100	100

Table 4: Data quality indicators for Zeolite A.