


Vertical STOE STADI MP



3 geometries in one diffractometer:

- Transmission/Debye-Scherrer
- Micro diffraction
- Bragg/Brentano

Movable tube housing, fixed goniometer

Monochromatic $K\alpha_1$ -radiation for all geometries

Mo, Cu and Co anodes

Various high and low temperature attachments

Stoe & Cie GmbH P.O.Box 101302 D-64213 Darmstadt

Phone: +49 6151 9887 0 Fax: +49 6151 9887 88

E-mail: stoe@stoe.com Homepage: <http://www.stoe.com>